



HaDeX - journey in the HDX-MS data


Weronika Puchała

About

- PhD student in Biophysics
- commercial background
- HDX-MS data analyst
- in the project from the beginning




Proteins!




But...


(a) Primary structure



(b) Secondary structure



(c) Tertiary structure




(d) Quaternary structure




HDX

H/D Exchange




Quench & Digest




©2013 David Weis

+ back-exchange!

MS



Largy, E. & Gabelica, V. Native Hydrogen/Deuterium Exchange Mass Spectrometry of Structured DNA Oligonucleotides. *bioRxiv* 848598 (2019) doi:10.1101/848598





The project and its difficulties


- new knowledge and context changes the workflow
- different fields of expertise
- terms changing the meaning
- local protocols
- advanced but narrow specialization



This Is Fine. Know Your Meme
<https://knowyourmeme.com/memes/this-is-fine>



Currently, in the project...

 1 Open
  0 Closed
 Sort ▾

HaDeX 1.5
 Updated 18 days ago


The aim of the project is to re-implement the background logic of the app to the package functions. General clean-up for easier and safer use of the package.
 ...

Latest release

 1.2
 ecdbd8a

Compare ▾


HaDeX 1.2


 **werpuc** released this on 1 Jun · 4 commits to master since this release


Changes from version 1.1 to 1.2:

- added tooltips to all of the plots
- improved the processing of peptides with modification
- renamed `kinetics` to `uptake curve`
- improved visibility of elements depending on the chosen calculation method

Assets 3

 [HaDeX_1.2.tar.gz](#)


 [Source code \(zip\)](#)


 [Source code \(tar.gz\)](#)


Feb 10, 2019 – Jul 20, 2020


Contributions: Commits ▾


Contributions to master, excluding merge commits





 **werpuc**
 247 commits 76,413 ++ 9,134 --




 **michbur**
 137 commits 28,313 ++ 45,475 --



 **DominikRafacz**
 16 commits 11,024 ++ 8,792 --




+



Our data

Protein	Start	End	Sequence	Modification	Fragment	Maxuptake	MHP	State	Exposure	File	z	RT	Inten	Center
1: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.000	KD_160527_CD160_sekw_05	1	3.232479	6592	1591.2584
2: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.000	KD_160527_CD160_sekw_05	2	3.238079	394066	796.3552
3: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.000	KD_160527_CD160_sekw_05	3	3.238759	173526	531.2633
4: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.001	KD_160527_CD160_IN_01	2	3.258598	232221	796.3634
5: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.001	KD_160527_CD160_IN_01	3	3.256844	110675	531.2849
6: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.167	KD_160530_CD160_10s_01	2	3.262774	99894	800.3610
7: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.167	KD_160530_CD160_10s_02	2	3.264558	117541	800.3852
8: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.167	KD_160530_CD160_10s_03	2	3.263795	90562	800.3682
9: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	0.167	KD_160530_CD160_10s_04	2	3.273656	66131	800.4242
10: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1.000	KD_160530_CD160_1min_01	2	3.262300	109301	800.7878
11: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1.000	KD_160530_CD160_1min_02	2	3.255346	112491	800.6889
12: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1.000	KD_160530_CD160_1min_03	2	3.254521	119940	800.8194
13: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1.000	KD_160530_CD160_1min_04	2	3.261150	105349	800.7712
14: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	5.000	KD_160530_CD160_5min_02	2	3.252980	139481	800.9808
15: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	5.000	KD_160530_CD160_5min_03	2	3.258085	108544	800.9552
16: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	5.000	KD_160530_CD160_5min_01	2	3.223357	329153	800.9635
17: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	5.000	KD_160530_CD160_5min_04	2	3.262770	68601	800.8628
18: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	25.000	KD_160531_CD160_25min_01	2	3.262472	116561	801.2475
19: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	25.000	KD_160531_CD160_25min_02	2	3.269511	62188	801.1659
20: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	25.000	KD_160531_CD160_25min_03	2	3.249441	95892	801.2836
21: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	25.000	KD_160531_CD160_25min_04	2	3.256427	92053	801.2446
22: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	120.000	KD_160531_CD160_2h_01	2	3.247047	140394	801.4678
23: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	120.000	KD_160531_CD160_2h_02	2	3.249840	115181	801.3489
24: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	120.000	KD_160531_CD160_2h_03	2	3.244796	147080	801.3217
25: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	120.000	KD_160531_CD160_2h_04	2	3.254010	143965	801.2839
26: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1440.000	KD_160601_CD160_24h_01	2	3.275502	96743	801.6348
27: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1440.000	KD_160601_CD160_24h_02	2	3.270558	100373	801.5983
28: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1440.000	KD_160601_CD160_24h_03	2	3.262421	82944	801.5474
29: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160	1440.000	KD_160601_CD160_24h_04	2	3.268210	86853	801.5782
30: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.000	KD_160527_CD160_sekw_05	1	3.232479	6592	1591.2584
31: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.000	KD_160527_CD160_sekw_05	2	3.238079	394066	796.3552
32: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.000	KD_160527_CD160_sekw_05	3	3.238759	173526	531.2633
33: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.001	KD_160527_CD160_IN_01	2	3.258598	232221	796.3634
34: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.001	KD_160527_CD160_IN_01	3	3.256844	110675	531.2849
35: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.167	KD_160531_HVEM_CD160_10s_01	2	3.251068	147694	800.2299
36: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.167	KD_160531_HVEM_CD160_10s_02	2	3.242326	178169	800.2436
37: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.167	KD_160531_HVEM_CD160_10s_03	2	3.254681	141814	800.2545
38: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	0.167	KD_160531_HVEM_CD160_10s_04	2	3.241787	194042	800.1728
39: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	1.000	KD_160531_HVEM_CD160_1min_01	2	3.250153	150273	800.6576
40: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	1.000	KD_160531_HVEM_CD160_1min_02	2	3.241649	190512	800.6585
41: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	1.000	KD_160531_HVEM_CD160_1min_03	2	3.249339	183445	800.4920
42: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	1.000	KD_160531_HVEM_CD160_1min_04	2	3.255107	100346	800.5765
43: db_CD160	1	15	INITSSASQEGTRLN	NA	NA	14	1590.8082	CD160_HVEM	5.000	KD_160531_HVEM_CD160_5min_01	2	3.254450	121446	800.6951

Welcome to HaDeXversum!



Welcome to HaDeX!

Thank you for using our tool.

Questions/feature requests/commercial applications: hadex@ibb.waw.pl

About

The HaDeX web server relies on functions from the **HaDeX R** package. For more information check our [documentation](#).

Authors

- Weronika Puchała
- Michał Burdukiewicz

Citation

Puchała W, Burdukiewicz M, Kistowski M, Dabrowska KA, Badaczewska-Dawid AE, Cysewski D and Dadlez M (2019). HaDeX: Analysis and Visualisation of Hydrogen/Deuterium Exchange Mass Spectrometry Data. R package version 1.0.

Funding

This work is supported by Foundation of Polish Science (TEAM TECH CORE FACILITY/2016-2/2 *Mass Spectrometry of Biopharmaceuticals - improved methodologies for qualitative, quantitative and structural characterization of drugs, proteinaceous drug targets and diagnostic molecules*).

Links

- Mass Spectrometry Laboratory, IBB PAS
- MIZ DataLab





Upload your file. Otherwise, you will see the example data.

Currently HaDeX is limited to `cluster` files from DynamX 3.0 or 2.0 and `tables` file from HDeXaminer

Accepted file extensions: .csv, .xls, .xlsx.

The supplied file should contain at least two repetitions of the measurement for the uncertainty to be calculated.

If the supplied file contains modified peptides, maximal exchange control cannot be applied.

Choose file:

Browse...

CK2_all results table.csv

Upload complete

File status:

Supplied file is valid. Detected data source: HDeXaminer. User action needed below!

Please be aware that loading data (including example file) may take a while. Be patient.

For the program to behave correctly, please make sure supplied file fulfills following requirements:

Show requirements

Settings

Values chosen here are propagated into all of the tables for coherent results.

Choose protein:

db_CD160

Choose maximal exchange control for chosen protein:

db_CD160 | CD160 | 1440

Choose D20 concentration [%]:

100

Move sequence start:

1

Correct sequence length:

132

Apply changes in `Input Data` tab.

File from HDeXaminer detected!

Some of the information from the data file requires your confirmation.

For the additional information on how the data from HDeXaminer is processed, check the requirements above. Keep in mind that the MHP value is generated based on the peptide sequence and therefore, may differ from actual value in case of the modifications.

FD timepoint [min]:

1440

Protein name:

Protein.State.1

States names:

Protein.State.1

Accepted confidence values:

☒ High

☒ Medium

☐ Low

Apply changes to continue

The calculated values of MPH might slightly differ based on data used and its precision.

Used amino mass data

Protein name


db_CD160

Reconstructed sequence

INITSSASQEGTRLNLICTVWHKKEEAEGFVVFLCKDRSGDCSPETSLKQLRLKRDPGIDGVGEISSQLMFTISQVTPLSHGTYYCCARSQKSGIRLQGHFFSILFxxxxxxxxxxxxxxxxxxFSHNEGTL

Name	Value
Length	132
Coverage	86.36%
Cys	5

- Hydro-
- ☒ Hydrophilic
 - ☒ Hydrophobic



[Save chart \(.svg\)](#)

Source: Kyte, J. and Doolittle, R.F. (1982). A simple method for displaying the hydropathic character of a protein. J Mol Biol 157, 105-132.

Choose state:

- ☒ CD160
- ☐ CD160_HVEM


Choose range:



Peptide Coverage

Data


Peptide coverage for db_CD160



Save chart (.svg)

Position Frequency

Data



Save chart (.svg)

Select parameters for the plot.

☐ Theoretical calculations

Choose values type:

- ☒ relative
☐ absolute

Comparison plot parameters:

Choose time parameters:

IN CHOSEN OUT

Choose states for comparison:

- ☒ CD160
☒ CD160_HVEM

Adjust colors

Woods plot parameters:

State 1 State 2

Confidence limit 1: Confidence limit 2:

Zoom:

Choose y range for comparison plot:



Choose y range for Woods plot:




Choose x range for both plots:



Adjust labels

Comparison plot


Data



Save chart (.svg)

Woods plot

Data



Save chart (.svg)

Select parameters for the plot.

☐ Theoretical calculations

Choose values type:

☒ relative

☐ absolute

Choose time parameters:

IN

☐ OUT

0.001

chosen control

Choose peptide:

Sequence	State	Start	End
INITSSASQEGTRLN	CD160	1	15
INITSSASQEGTRLN	CD160_HVEM	1	15
LICTWW	CD160	16	21
LICTWW	CD160_HVEM	16	21
LICTWWHKKEEAEG	CD160	16	29
LICTWWHKKEEAEG	CD160_HVEM	16	29
ICTVWHKKEEAEG	CD160	17	29
ICTVWHKKEEAEG	CD160_HVEM	17	29
PVVFL	CD160	30	34
PVVFL	CD160_HVEM	30	34

All

All

All

All

Showing 1 to 10 of 82 entries

1

2

3

4

5

...

9

Next

Reset chosen peptides

Choose y range for the plot:




Adjust labels

Uptake curve

Data

Uptake curve for chosen peptides for db_CD160



Save chart (.svg)

Select parameters for the plot.

The plot can be rendered only for relative experimental values.

Choose in time:

0.001

Choose time:

0.167

State 1


CD160

State 2

CD160_HVEM

Quality control plot

Data



Save chart (.svg)

This function plots the change in the uncertainty of deuteration levels as a function of incubation time. The uncertainty is averaged over all peptides available at a given time point in a selected state. This chart has a double function: firstly, it allows checking if the measurement uncertainty is decreasing over time (which is the expected behavior). Secondly, it helps to plan the appropriate incubation length for the tested protein (whether we obtain the desired data reliability values).

Name	Value
HDX time course	8
Number of peptides	41
Sequence coverage	86.36%
Average peptide length	14.6829
Redundancy	4.5606
Replicates	4
Significant differences in HDX	0.98 0.99
All	All

[Excel](#)
[PDF](#)

Showing 1 to 7 of 7 entries

This table was adapted from the HDX Data Summary (Masson, G.R., Burke, J.E., Ahn, N.G., Anand, G.S., Borchers, C., Brier, S., Bou-Assaf, G.M., Engen, J.R., Englander, S.W., Faber, J., et al. (2019), Recommendations for performing, interpreting and reporting hydrogen deuterium exchange mass spectrometry (HDX-MS) experiments. Nat Methods 16, 595–602).

Choose items for report:

- | | |
|---|---|
| <input checked="" type="checkbox"/> Position Frequency | <input type="checkbox"/> Position Frequency Data |
| <input checked="" type="checkbox"/> Peptide Coverage | <input type="checkbox"/> Peptide Coverage Data |
| <input checked="" type="checkbox"/> Comparison Plot | <input type="checkbox"/> Comparison Plot Data |
| <input checked="" type="checkbox"/> Theoretical Comparison Plot | <input type="checkbox"/> Theoretical Comparison Plot Data |
| <input checked="" type="checkbox"/> Woods Plot | <input type="checkbox"/> Woods Plot Data |
| <input checked="" type="checkbox"/> Theoretical Woods Plot | <input type="checkbox"/> Theoretical Woods Plot Data |
| <input type="checkbox"/> Uptake Curve | <input type="checkbox"/> Uptake Curve Data |
| <input type="checkbox"/> Theoretical Uptake Curve | <input type="checkbox"/> Theoretical Uptake Curve Data |
| <input type="checkbox"/> Quality Control Plot | <input type="checkbox"/> Quality Control Plot Data |

Elements chosen for report have the same parameters as chosen in panel e.g. axis range and title. Adjust parameters for plots as needed in the report.

[Create report!](#)

Keeping up with the competitors



MEMHDX

Finally!

ACCEPTED MANUSCRIPT

HaDeX: an R package and web-server for analysis of data from hydrogen-deuterium exchange mass spectrometry experiments

Weronika Puchała, Michał Burdukiewicz, Michał Kistowski,
Katarzyna A Dkšbrowska, Aleksandra E Badaczewska-Dawid,
Dominik Cysewski ✉, Michał Dadlez ✉ Author Notes

Bioinformatics, btaa587, <https://doi.org/10.1093/bioinformatics/btaa587>

Published: 24 June 2020 **Article history** ▼

Thanks!